



**Thermo**  
SCIENTIFIC

It's time to  
be sure



Thermo Scientific  
SureTect Real-Time  
PCR System



# It's better to be sure

Now, you can perform every food safety test confidently with the **Thermo Scientific™ SureTect™ Real-Time PCR System**. Built on proven PCR technology and backed by world-class service and support, the SureTect™ System is designed to quickly and accurately detect microorganisms in a broad range of foods and associated samples. This unique solution combines speed and performance in an easy-to-use, cost-effective platform—giving you results you can be sure about.



**Thermo**  
SCIENTIFIC



# It's easy to be sure

At Thermo Fisher Scientific, we understand microbiologists. We know that you need fast, reliable molecular technology to get through the day-to-day challenges. That's why we continue to invest in bringing you complete solutions, like the SureTest Real-Time PCR System. If you want high quality results without the stress, all it takes is three simple steps.

1.



## Prepare

Reduce hands-on time and eliminate variance with pre-filled assay tubes, and fast, simple lysis in less than 20 minutes

## Run

Up to five 24-well PikoReal instruments can be run from the same PC, to accommodate laboratories of all sizes

2.



3.



## Read

Intuitive software enables rapid, straightforward interpretation and traceability of results



## Why choose SureTect?

### Single enrichment step

no secondary enrichment or regrowth required, for faster, simpler testing

### Pre-filled lysis tubes

for maximum convenience and reliable, consistent cell lysis

### Tableted PCR reagents

eliminates PCR reagent pipetting, to enhance reproducibility and minimize hands-on time

### Streamlined protocol

facilitates processing of multiple assays in the same run, for maximum efficiency

### Simple, intuitive software

for straightforward training, quick set-up, and simple tracking of results

### Automatic data interpretation

reliably reports test results as positive or negative

## Unmatched service & support

Thermo Fisher Scientific experts have been supporting the global microbiology community for over 150 years. Our service and support teams are comprised of highly-experienced scientists who are devoted to offering best-in-class technical support.

Our focus on providing quality products, on-time delivery and superior support is matched by our commitment to provide complete solutions that meet your microbiological testing needs.

For more information, please contact your local Thermo Scientific representative by visiting [www.thermoscientific.com/SureTect](http://www.thermoscientific.com/SureTect).

[thermoscientific.com/SureTect](http://thermoscientific.com/SureTect)

© 2013 Thermo Fisher Scientific Inc. All rights reserved.  
All trademarks are property of Thermo Fisher Scientific Inc., and its subsidiaries.

**Thermo**  
SCIENTIFIC  
Part of Thermo Fisher Scientific