

Technical Bulletin

Xpect[®] Flu A&B IFU Update Effective August 18, 2009: Additional transport media validated for use with Xpect[®] Flu A&B.

Xpect[®] Flu A&B (R24600)

The following transport media have been evaluated and found to be compatible with the Xpect[®] Flu A&B Test Kit.

- Amies Medium
 - Bartels Viral Transport Medium
 - BD™ Universal Viral Transport System (**updated**)
 - Cary Blair Medium
 - Copan Universal Transport Medium (**updated**)
 - Earle's Minimum Essential Medium (EMEM)
 - EMEM with 1% Bovine Serum Albumin
 - EMEM with 1% Lactalbumin Hydrolysate
 - Hank's Balanced Salt Solution
 - Liquid Stuarts Medium
 - M4[®]
 - M4RT[®]
 - M5[®]
 - M6™ (**updated**)
 - Phosphate Buffered Saline (PBS)
 - PBS with 0.5% Bovine Serum Albumin
 - PBS with 0.5% Gelatin
 - Saline (normal)
 - Sucrose Phosphate
 - Tryptic Soy Broth with 0.5% Bovine Serum Albumin
 - Tryptic Soy Broth with 0.5% Gelatin
 - Veal Infusion Broth
 - Veal Infusion Broth with 0.5% Bovine Serum Albumin
-

Please contact Remel Technical Service at 800.255.6730 or itechnicalsupport@remel.com if you have any questions or require any assistance with Remel products.