



## More antifungal options lower cost per test

Thermo Scientific Sensititre YeastOne plates, including the only FDA-cleared broth microdilution plate for antifungal susceptibility testing

- Colorimetric Thermo Scientific™ alamarBlue™ agent provides reliable, easy and consistent endpoint determination with visual read options
- Dual isolate, multi-antifungal IVD format reduces cost per test for *Candida* spp.
- Single isolate, nine antifungal, research-use-only format includes anidulafungin\* and micafungin\* for yeast and filamentous fungi
- Room-temperature storage eliminates inventory control concerns
- Individual packaging allows laboratory to test one plate at a time with no waste
- Inclusive on-scale QC ranges provide immediate quality assurance of testing methodology
- 24-hour incubation time enables timely results for clinicians

\*For research use only in the U.S. Not for use in diagnostic procedures.

**Thermo**  
SCIENTIFIC

## Sensititre YeastOne Plates:


Antimicrobial:	Part #: YO2IVD (IVD) Dilution Range (µg/mL):	YO9*,** (RUO)	YO10***
Amphotericin B		0.12-8	0.12-8
5-Flucytosine	0.03-32	0.06-64	0.06-64
Anidulafungin		0.015-8	0.015-8
Caspofungin	0.015-8	0.008-8	0.008-8
Micafungin		0.008-8	0.008-8
Fluconazole	0.25-64	0.12-256	0.12-256
Itraconazole	0.03-2	0.015-16	0.015-16
Posaconazole		0.008-8	0.008-8
Voriconazole	0.008-4	0.008-8	0.008-8

### Related products:

Sensititre YeastOne IVD Plate  
 Sensititre YeastOne Plate w/ Anidulafungin and Micafungin  
 Sensititre YeastOne CE Plate w/ Anidulafungin and Micafungin  
 Sensititre YeastOne Broth, 11mL; 10/box  
 0.5 McFarland Standard  
 Demineralized Water, 5mL; 100/box  
 Susceptibility Organism Kit, Sensititre YeastOne Training Set; 10 organisms  
 Thermo Scientific Culti-Loops QC Organism Set for Sensititre YeastOne plates

### Part #:

YO2IVD  
 YO9\*,\*\*  
 YO10\*\*\*  
 Y3462  
 E1041  
 T3339  
 T5102\*,\*\*  
 R4687027



Sensititre YeastOne plate image on the Thermo Scientific Sensititre Vizion Digital MIC Viewing System.

\*For research use only in the U.S. Not for use in diagnostic procedures.

\*\*Available in the U.S. only.

\*\*\*Not available in the U.S.

© 2012 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**USA and Canada** +1 800 871 8909  
**All Other Inquiries** +44 1256 841144

**Technical Support**  
**USA and Canada** +1 800 642 7029  
**All Other Inquiries** +44 1342 318777

**Thermo**  
 S C I E N T I F I C

Part of Thermo Fisher Scientific

LT2024A/CS/02/12