



RapID™ SS/u

Report Form

Reference # / No. de référence / Referenz-Nr. / Riferimento N. / N° de referencia _____

Date / Date / Datum / Data / Fecha _____

Tech / Tech / Techn. / Tech / Tec _____

Source / Source / Quelle / Origine / Origen _____

Oxidase / Oxydase / Oxidase / Ossidasi / Oxidasa _____

Growth Media / Milieu de croissance / Nährböden / Terreno di crescita / Medios de crecimiento _____

Reagent/ Réactif/ Reagenz/ Reagente/ Reactivo	None / Aucun / Keine / Nessuno / Ninguno						IDS RapID SS/u Reagent / Réactif RapID SS/u / RapID SS/u Reagens /RapID SS/u Reagent / Reactivo RapID SS/u			Spot indole					
Positive Reactions Réactions positives Positive Reaktionen Reazioni positive Reacciones positivas	Medium or bright yellow Jaune moyen ou vif Mittleres oder starkes Gelb Giallo di intensità media o brillante Amarillo medio o brillante						Dark red, red, red-orange or orange Rouge foncé, rouge, rouge orangé ou orange Dunkelrot, Rot, Rot-Orange oder orange Rosso scuro, rosso, rosso-arancio o arancio Rojo oscuro, rojo, rojo-naranja o naranja			Purple, violet, red, or dark pink Violacé, violet, rouge ou rose soutenu Purpur, Violett, Rot oder Dunkelrosa Porpora, violetto, rosso o rosa scuro Morado, violeta, rojo o rosa oscuro			Black, brown, or tan Noir, marron ou légère coloration Schwarz, Braun oder dunkel getönt Nero, marrone o beige Negro, marrón o tostado		
Cavity # / No. cavité / Kammer-Nr. / Cavià N. / N° de cavidad	1	2	3	4	5	6	7	8	9	10	7				
Test Code / Code du test / Testcode / Codice esame Código de prueba	GMS	ONPG	G1	G2	G3	PHS	URE	A1	A2	A3	IND	OXI			
Value / Valeur / Wert Valore / Valor	1	2	4	1	2	4	1	2	4	1	2	4			
Result / Résultat Ergebnis / Risultato Resultado															
Value Total Total des valeurs Gesamtwert Valore totale Valor total															

IDENTIFICATION / IDENTIFICATION / IDENTIFIZIERUNG / IDENTIFICAZIONE /
IDENTIFICACIÓN _____

